


# ComVi Pyloric/Duodenal Stent [Flare-Type]

for intrinsic and/or extrinsic malignant pyloric/duodenal stricture



## Features

- **Triple layered construction:** A biocompatible PTFE membrane tube is held between the inner and outer mesh
  - An uncovered proximal flare end to reduce migration
  - The PTFE membrane prevents the risk of tissue invasions
- **Large diameters (up to 26mm)** are loaded into a 10.5Fr delivery system
- **Radiopaque marker:** 3 (three) at the Flared end and 2 (two) on both ends of the body

# Ordering Information

OTW (Over The Wire)				
Both Ends Bare Type				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PC1806BP	18	6	14	135
PC1808BP		8		
PC1810BP		10		
PC1812BP		12		
PC2006BP	20	6	14	135
PC2008BP		8		
PC2010BP		10		
PC2012BP		12		
PC2206BP	22	6	16	135
PC2208BP		8		
PC2210BP		10		
PC2212BP		12		
PC2406BP	24	6	16	135
PC2408BP		8		
PC2410BP		10		
PC2412BP		12		
PC2606BP	26	6	18	135
PC2608BP		8		
PC2610BP		10		
PC2612BP		12		
PC2806BP	28	6	18	135
PC2808BP		8		
PC2810BP		10		
PC2812BP		12		

TTS (Through The Scope)				
Distal End Bare Type				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PCT1806P	18	6	10.5	180
PCT1808P		8		
PCT1810P		10		
PCT1812P		12		
PCT2006P	20	6	10.5	180
PCT2008P		8		
PCT2010P		10		
PCT2012P		12		
PCT2206P	22	6	10.5	180
PCT2208P		8		
PCT2210P		10		
PCT2212P		12		
PCT2406P	24	6	10.5	180
PCT2408P		8		
PCT2410P		10		
PCT2412P		12		
PCT2606P	26	6	10.5	180
PCT2608P		8		
PCT2610P		10		
PCT2612P		12		

# Ordering Information

TTS (Through The Scope)				
Distal End Bare Type				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PCT1806P-12	18	6	10.5	120
PCT1808P-12		8		
PCT1810P-12		10		
PCT1812P-12		12		
PCT2006P-12	20	6	10.5	120
PCT2008P-12		8		
PCT2010P-12		10		
PCT2012P-12		12		
PCT2206P-12	22	6	10.5	120
PCT2208P-12		8		
PCT2210P-12		10		
PCT2212P-12		12		
PCT2406P-12	24	6	10.5	120
PCT2408P-12		8		
PCT2410P-12		10		
PCT2412P-12		12		
PCT2606P-12	26	6	10.5	120
PCT2608P-12		8		
PCT2610P-12		10		
PCT2612P-12		12		

TTS (Through The Scope)				
Distal End Bare Type				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PCT1806P-22	18	6	10.5	220
PCT1808P-22		8		
PCT1810P-22		10		
PCT1812P-22		12		
PCT2006P-22	20	6	10.5	220
PCT2008P-22		8		
PCT2010P-22		10		
PCT2012P-22		12		
PCT2206P-22	22	6	10.5	220
PCT2208P-22		8		
PCT2210P-22		10		
PCT2212P-22		12		
PCT2406P-22	24	6	10.5	220
PCT2408P-22		8		
PCT2410P-22		10		
PCT2412P-22		12		
PCT2606P-22	26	6	10.5	220
PCT2608P-22		8		
PCT2610P-22		10		
PCT2612P-22		12		